



IMMUNO REAGENTS

Encephalitozoon cuniculi whole cell antigen suspension



Features

- Whole cell suspension
- High antigenicity

Product description

The *Encephalitozoon cuniculi* (*E. cuniculi*) antigen suspension consists of about 3×10^7 spores per ml. The spores, killed by heating, are washed and suspended in phosphate buffered saline solution (PBS) containing 0.1% formalin. The *E. cuniculi* antigen suspension is to be used together with carbon suspension in Medicago's CIA test.

For laboratory use only.

Applications

- Diagnostics
- Immunology

Specifications

Formulation	1x PBS buffer, 0.1% formaldehyde
Concentration	3×10^7 spores per ml
Storage	2°C to 8°C
Shelf life	Three years after production date

Shipping and storage

The product is shipped at room temperature and antigenicity will last for three years or more. Store at 2°C to 8°C.

Certifications

Medicago's laboratories and manufacturing site in Uppsala are ISO 9001:2008 and ISO 13485:2003 certified. Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



Ordering information

Article no.	Product name	Pack size
18-2002-3	<i>E. cuniculi</i> whole cell antigen suspension	3 ml