



IMMUNO REAGENTS

Rabbit antiserum against *Toxoplasma gondii*



Features

- A rabbit serum that can be used as a positive control in CIA procedures

Product description

Toxoplasma gondii is a species of parasitic protozoa. The primary host of *T. gondii* is the cat, but the parasite can be carried by all known mammals. *T. gondii* is the causative agent of Toxoplasmosis, a disease that is usually minor and self-limiting but that can have serious or even fatal effects on a fetus whose mother first contracts the disease during pregnancy, or on an immuno-compromised human or cat.

For laboratory use only.

Applications

- Veterinary procedures
- A rabbit serum that can be used as a positive control in immunoassays

Directions for use

Dilute serum 1:20 or more in PBS buffer. Use lab coat and gloves when handling. Avoid repeated freeze-thawing.

Shipping and storage

The product is shipped at 2-8°C, store serum below -18°C.

Certifications

Medicago's laboratories and manufacturing site in Uppsala are ISO 9001:2008 and ISO 13485:2003 certified. Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



Specifications

Content	Serum from rabbits, immunized with <i>T.gondii</i> antigen
Storage	Below -18°C
Shelf life	12 months undiluted

Ordering information

Article no.	Product name	Pack size
18-5001-100	Rabbit antiserum against <i>T.gondii</i>	0,1 ml
18-5001-1	Rabbit antiserum against <i>T.gondii</i>	1 ml